



PERSONAL PROTECTIVE EQUIPMENT (PPE)

FACE SHIELD – Model: AEPPE 2304

AEON ELECTRONICS PVT LTD
Complete Manufacturing Solution Provider

Regd. Office: No 25, 2nd Stage, 6th Main Road,
J C Industrial Layout, Kanakapura Road,
Bangalore – 560 062, India
www.aeontele.com

Suresh Babu: +91 9343780336
vasanth@aeontele.com
Kumar@aeontele.com



MODEL: AEPPE 2304

The Product

- Made of 350 microns clear transparent sheet
- Can be used with safety Glass and Mask
- Adjustable Velcro tape to fit any head size
- Forehead cushion for comfort and to accommodate wearing safety glasses
- Comfortable for long usage
- Sterilized packing

Application

- Used for daily protection – Both indoor and Outdoor
Examples include hospitals, retail shop, police check post, petrol pumps, bank, factory, office, etc.

User Friendly Product

- Unisex and Unisize head clip that's lightweight
- User Friendly Product as it is very easy to wear

Technical Details

- Sheet Thickness: 360 GSM
- Number of units per pack: 60
- Size of transparent Sheet: 300 mm x 240 mm
- Size of Velcro band: 750 mm

Endurance Parameters

- Impact Resistance: Good
- Scratch Resistance: Good
- Heat Resistance: Good
- Optics: Clear and Good



Forehead Cushion foam

Velcro Band for head fitment

Transparent sheet

Instructions for Use

- Before handling the package, wash or sanitize your hands
- Remove the Face Shield from the sealed cover carefully
- Do not keep the face shield on any surface after removing from the sealed cover
- Place the forehead pad and use Velcro band to comfortably fit to your head

Note: Depending on the usage, Face mask and safety goggles should be worn before using Face Shield

Caution

This face shield is developed understanding the importance of comfort, flexibility, fit and ease of use. These protects eyes, nose and mouth from splash of droplets. Face shield are usually used with conjunction with other personal protective equipment such as face mask, goggles based on the hazard potential. Please take necessary advice depending on your usage circumstances.

Pls Dispose like any other Medical plastic Hazard material and Sanitize or wash your hands