



## **PERSONAL PROTECTIVE EQUIPMENT (PPE)**

---

**FACE SHIELD – Model: AEPPE 2404**

**AEON ELECTRONICS PVT LTD**  
*Complete Manufacturing Solution Provider*

**Regd. Office:** No 25, 2nd Stage, 6th Main Road,  
J C Industrial Layout, Kanakapura Road,  
Bangalore – 560 062, India  
[www.aeontele.com](http://www.aeontele.com)

Suresh Babu: +91 9343780336  
vasanth@aeontele.com  
Kumar@aeontele.com



MODEL: AEPPE 2404

### The Product

- Each FACE SHIELD includes a reusable head clip and one transparent sheet
- The clear transparent sheet is made from high quality PET
- The reusable head clips are made from high quality plastic raw material

### Application

- Used for daily protection – Both indoor and Outdoor  
Examples include hospitals, retail shop, police check post, petrol pumps, bank, factory, office, etc.

### User Friendly Product

- Unisex and Unisize head clip that's lightweight
- Flexible and reusable multiple times with sanitization
- Additional transparent sheets are readily available and can be purchased separately
- The design is convenient for the users to replace the transparent sheet at regular intervals

### Technical Details

- Standard: ANSI Z87.1-2010 (for splash of liquid droplets)
- Weight per face shield: 35 grams
- Number of units per pack: 100
- Size of transparent Sheet: 340 mm x 220 mm

### Endurance Parameters

- Impact Resistance: Good
- Scratch Resistance: Good
- Heat Resistance: Good
- Optics: Clear and Good





#### Face Shield

#### ASSEMBLY INSTRUCTIONS

#### SAFETY INSTRUCTIONS:

1. Before handling the package, wash your hands with sanitizer.
2. Clean the head band with sanitizer and use cleaned area for assembly



### Instructions for Reuse

- Remove the Transparent Sheet carefully
- Sanitize the head clip with soap water or suitable sanitizer and dry it immediately
- Replace the transparent film and it's ready to use

### Caution

This face shield is developed understanding the importance of comfort, flexibility, fit and ease of use. These protects eyes, nose and mouth from splash of droplets. Face shield are usually used with conjunction with other personal protective equipment such as face mask, goggles based on the hazard potential. Please take necessary advice depending on your usage circumstances.

Pls Dispose like any other Medical plastic Hazard material and Sanitize or wash your hands